

co₂ntrol[®]


Sofnolime



Sofnolime[®] pre-filled canister for closed anaesthesia

Applications

co₂ntrol is a pre-filled disposable CO₂ absorber designed for use with Dräger anaesthesia systems. Used for the removal of CO₂ in breathing systems where the pressure-drop must be negligible, the main component is calcium hydroxide. Contains a white to violet indicator to show absorption of CO₂. For use with Dräger workstations fitted with the Dräger CLIC[®] adapter MX50090.



Packing and transport

co₂ntrol is available filled with Sofnolime (1.15kg) and Sofnolime SoLo (1.0kg) medical USP grade soda lime.

co₂ntrol canisters are packed 6 in a box; 27 boxes per pallet (standard). The product has a shelf-life of 2 years.

Sofnolime specification

Particle size distribution¹

Greater than 8 mm	%	nil	nil
Between 2 and 5 mm	%	78.0	Min
Up to 0.425 mm	%	0.5	Max

Friability ² (equivalent to USP hardness)	%	95.0	Typical
Moisture absorption	%	7.5	Max
USP CO ₂ activity ³	%	19.0	Min

Properties

Water content	%	12 - 19
---------------	---	---------

CO₂ performance (canister)

130 litres of carbon dioxide per canister at NTP (20°C and 1 bar @ 0.475 litre CO₂ per min until 0.5% CO₂ breakthrough)

Quality

Molecular Products Ltd's aim is to manufacture chemical products which completely satisfy the needs of our customers. All products are rigorously tested to ensure conformance to the specification.

Our activities comply to the requirements of ISO 9001 and ISO 13485. Soda lime, as a class IIa medical device, is manufactured under ISO 13485 procedures.

Notes

1. Particle size measured by optical method
2. Friability measured by ball mill and optical method
3. Activity is determined by passing carbon dioxide through a 'U' tube containing approximately 10g of Sofnolime at a rate of 75ml per minute for a period of 20 minutes. The activity is the weight gain, expressed as a percentage
4. Particle size and friability methods are validated against USP test method
5. Particle size distribution is controlled to provide standardised pressure drop
6. Determined using internal test method TM153. Test conditions: 500ml tidal volume, breathing rate of 20 breaths per min (bpm) and an I:E ratio to 1:2

Molecular Products Limited

Parkway, Harlow Business Park
Harlow, Essex, CM19 5FR, UK

T +44 (0)1279 445111
F +44 (0)1279 401231

E sales@molprod.com
W www.molecularproducts.com

