# SODASORB®

## CO<sub>2</sub> ABSORBENT

### STANDARD SODASORB

1. 使用说明 (CHIN)
2. NÁVOD K POUŽITÍ (CZ)
3. VERWENDUNGSHINWEISE (DE)
4. INSTRUKTIONER FOR BRUG (DK)
5. INSTRUCCIONES PRELIMINARES PARA EL USO (ES)
6. KÄYTTÖOHJEET (FIN)
7. INSTRUCTIONS D'UTILISATION (FR)
8. INSTRUCTIONS FOR USE (GB)
9. ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ (GR)
10. প্রযোগ নির্দেশ (HIN)
11. HASZNÁLATI UTASÍTÁSOK (HU)
12. ISTRUZIONI PER L'USO (IT)
13. 使用説明書 (J)
14. BRUKSANVISNING (N)
15. GEBRUIKSAANWIJZINGEN (NL)
16. MODO DE EMPREGO (P)
17. INSTRUKCJA UŻYCIA (PL)
18. Инструкция по применению (RUS)
19. BRUKSANVISNING (SW)
20. KULLANIM TALİMATLARI (TR)
21. DISCLAIMER

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# SODASORB®

## Carbon Dioxide Absorbent

**HENKEL Corporation aka  
HENKEL Technology Corporation**  
6050 West 51st Street,  
Chicago, IL, 60638, United States  
Tel : +1 200-GET-SODA  
Tel : +1 800-438-7632

### **GENERAL DESCRIPTION & PACKAGING**

SODASORB® absorbent is intended for medical purposes to remove carbon dioxide from gases in the breathing circuit in systems such as anaesthesia and respiratory therapy equipment.

**Caution: Federal law restricts this device to sale by or on the order of a physician (U.S.A.)**

SODASORB® absorbent is available in a variety of package configurations:

Plastic Jug (4.5 kg), carton of 5 pieces. Shelf life = 3 years  
Plastic Jug (20 kg). Shelf life = 3 years

### **PREPARATION**

Opening / unwrapping:

- Carefully open containers to avoid any spills.

Filling and changing:

- Anesthesia equipment absorbent canisters should be cleaned when the absorbent is changed, with particular attention paid to any screens, as they are susceptible to obstruction.
- Absorbents should always be handled gently to avoid fragmentation and dust formation. When a canister is emptied, care should be taken to remove dust particles,
- The canister should always be filled with care. Loose-fill absorbent should be poured slowly into the canister while the canister is rotated, stopping occasionally to tap the sides. The canister should be filled completely but not overfilled. A small space should be left at the top to promote even gas flow through the canister. The upper layer of absorbent should be level.

Pre-check:

- It is the responsibility of operating room personnel to perform any necessary equipment pre-checks.

### **USE**

- SODASORB® absorbent is for use with oxygen, nitrous oxide, halothane, enflurane, desflurane, sevoflurane and isoflurane. Before using any other anaesthetics agent, obtain approval from the anaesthetic manufacturer

**Caution: do not use SODASORB® absorbent with chloroform or trichloroethylene**

- Before use, check the breathing circuit for gas flow restrictions.
- **SODASORB® absorbent contains an indicator that turns from white to purple as the absorption capacity approaches exhaustion.** The absorption of carbon dioxide is indicated by a gradual color change in the direction of air flow through the absorbent. The intensity of color change may vary from one

application to another. Color change is an adjunct guide, the anaesthetist should rely primarily on CO<sub>2</sub> monitors, as well as time and volume calculations, to determine absorptive capacity.

**Caution: color change is not permanent, Discard product as soon as color change has occurred.**

- Operating room personnel shall determine how much absorptive capacity is left in the SODASORB® absorbent immediately at the end of the surgical procedure.
- Oxygen flush duration should not be overused. Too much high flushing of oxygen induces dehydration of SODASORB® absorbent and inhibits its ability to absorb carbon dioxide.
- SODASORB® absorbent is not an anti-microbial filter or inhibitor. It is the responsibility of the surgeon to take appropriate measures regarding patient cross contamination.

**Caution: there is potential for some portion of SODASORB® absorbent dust to migrate into the breathing circuit; the possibility of occurrence varies based on the conditions of use, handling, application, and type of anesthesia machine.**

- Filters are recommended to protect the patient from contaminants.
- SODASORB® absorbent should be used within its shelf life. The expiration date is printed on the packaging unit.
- Re-sealable containers, once opened, must be properly resealed to maintain product integrity and should in all cases be used within one month of opening.

### **STORAGE & DISPOSING**

Storage and Handling:

**Caution: irritant**

- Avoid contact with eyes, skin and clothing.
- Store product in an environment that minimizes prolonged exposure to UV light.
- Avoid creating airborne dust, especially when pouring or handling loose material.
- Keep from freezing, which may reduce CO<sub>2</sub> absorption performance and increase fines.
- Handle with care, especially when placing product in a canister, to avoid breathing dust, which can cause irritation.
- Consult MSDS for additional information.

Disposal:

- Disposal of used SODASORB® absorbent must be performed by the end user in compliance with all applicable laws and local regulations. Specific disposal methods may vary with the chemical (e.g., anaesthetic agents, acid-base vapours, etc.) or biological agents used with the absorbent by the end user.
- Keep all SODASORB® absorbent that has been used with highly flammable anaesthetic away from heat, sparks and open flames, as residual amounts of these materials will be present.

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**GB**

The information presented herein is derived from our testing and experience. It is offered, free of charge, for your consideration, investigation and verification. Since operating conditions vary significantly, and since they are not under our control, we disclaim any and all warranties on the results which might be obtained from the use of our products. You should make no assumption that all safety or environmental protection measures are indicated or that other measures may not be required.

SODASORB®, SODASORB® LF are trademarks, registered in the United States and/or other countries, of HENKEL Corporation. This trademark list has been compiled using available published information as of the publication date of this document and may not accurately reflect current trademark ownership or status. © Copyright 2017 HENKEL Corporation. All rights reserved.

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SODASORB®, SODASORB® LF 是 HENKEL Corporation 公司在美国和/或其他国家的注册商标。本商标名单采用本文件发布日起已公布的信息进行编撰，故可能无法准确反映当前商标的所有权或状态。版权所有© Copyright 2017 HENKEL Corporation. All rights reserved.
SODASORB®, SODASORB® LF jsou obchodní známky společnosti HENKEL Corporation, registrované ve Spojených státech a/nebo jiných zemích. Tento seznam obchodních známek byl vytvořen s pomocí dostupných zveřejněných informací k datu zveřejnění tohoto dokumentu, nemusí přesně odrážet aktuální vlas-tnický nebo statut obchodních známek. © Copyright 2017 HENKEL Corporation. Všechna práva vyhrazena.
SODASORB®, SODASORB® LF er varemærker registreret i USA og/eller andre lande for HENKEL Corporation. Denne varemærkeliste er udfulgt ved anvendelse af til rádgående udgivet information fra og med udgivelsesdatoen af dette dokument og er muligvis ikke en præcis status over ejerskabet til varemærket. © Copyright 2017 HENKEL Corporation. Alle rettigheder forbeholdt.
SODASORB®, SODASORB® LF son marcas comerciales de HENKEL Corporation registradas en los Estados Unidos y otros países. Esta lista de marcas registradas se ha recopilado utilizando información publicada a partir de la fecha de publicación de este documento y puede no ser un reflejo fiel de la propiedad y estado de las marcas. © Copyright 2017 HENKEL Corporation. Reservados todos los derechos.
SODASORB®, SODASORB® LF ovat tuotemerkejä, jotka HENKEL Corporation on rekisteröinyt Yhdysvalloissa ja/tai muissa maissa. Tämä tuotemerkkilista on laadittu käytännällä julkisia tietoja lukien tämän asiakirjan julkaispäivästä ja se ei välittämättä kuvata oikein tuotemerkin nykyisistä omistajuutta tai statusta. © Copyright 2017 HENKEL Corporation. Kaikki oikeudet pidätetään.
SODASORB®, SODASORB® LF a HENKEL Corporation Egyesült Államokban és/vagy más országokban bejegyzett védjegyei. A védjegyek listája a dokumentum műegyenlőségi döntéssel időpontjában nyilvánosságra hozott információk alapul, és elközelhető, hogy nem felel meg az aktuális védjegy státusz zoknak. © Copyright 2017 HENKEL Corporation. Minden jog fenntartva.