

RUGGED OXYGEN GENERATOR (ROG)

TECHNICAL PAPER

By Andrew Maier



ABSTRACT

The capability to deliver medical oxygen during prehospital combat care remains a critical gap in Large-Scale Combat Operations (LSCO) and other degraded environments where power, logistics, and safety constraints limit the use of conventional oxygen systems.

Compressed gas cylinders present fragmentation and explosive hazards under kinetic threat, require complex sustainment, and are often withheld from the most dangerous terrain. Portable oxygen concentrators (POCs), while effective in clinical settings, are dependent on electrical power and experience reductions in flow and oxygen purity at altitude, under thermal stress, and during high-demand ventilation. The absence of a safe, self-contained oxygen capability diminishes survivability for casualties requiring airway and ventilatory support during the most time-sensitive phase of care.

The Rugged Oxygen Generator (ROG) addresses this gap through a solid-state chemical oxygen generation system capable of delivering a minimum of 96% oxygen at ≥ 6 L/min for a minimum of 15 minutes, without external power, pressurization, or maintenance. The ROG has completed MIL-STD-810G environmental tests, is FDA 510(k) cleared, and UN3356 certified for global transport.

This capability directly supports Tactical Combat Casualty Care (TCCC) guidelines and enhances survivability when treatment is needed at point of injury due to enemy threat, distance, or airspace denial.

TABLE OF CONTENTS

Abstract.....	i
1. Introduction / Background.....	3
2. Problem Statement / Capability Gap.....	3
3. Operational Impact of the Gap.....	4
Figure 1. Hypoxic Deterioration Sequence in Combat Trauma.....	5
4. Operational & Physiological Consequences of Hypoxia.....	5
5. Evaluation Criteria for an Effective Solution.....	6
6. Solution Overview: ROG Rugged Oxygen Generator.....	7
7. Testing & Validation.....	9
Table 1. Summary of MIL-STD-810G Tests Performed.....	9
8. Comparative Analysis.....	10
Table 2. Capability Comparison of Oxygen Systems.....	11
9. Cross-Domain Medical Applications.....	11
10. Acquisition & Sustainment Pathway.....	12
11. References.....	14

1. INTRODUCTION / BACKGROUND

Oxygen is a primary therapeutic intervention in combat trauma care, particularly for casualties with traumatic brain injury (TBI), airway compromise, burns, blast lung, and thoracic injury. Hypoxia accelerates secondary brain injury and multi-organ failure, and even brief interruptions in oxygen delivery increase mortality and long-term neurological morbidity. Short hypoxic episodes during the prehospital phase are independently associated with poorer functional outcomes and higher mortality (Spaite et al., 2017).

During U.S. combat operations in Iraq and Afghanistan, approximately 24.3 percent of battlefield deaths were classified as potentially survivable, and 87 percent of those occurred before arrival at a medical treatment facility (Eastridge et al., 2012). The leading causes were hemorrhage and airway-respiratory failure - conditions for which early oxygenation is essential to prevent irreversible physiological deterioration.

Future large-scale combat operations (LSCO) are expected to impose extended evacuation timelines because of distributed maneuver, austere terrain, and contested airspace. Medical forces must anticipate evacuation delays that exceed doctrinal goals, increasing dependence on prolonged prehospital care. Current oxygen systems are ill-suited to these environments: pressurized oxygen cylinders are withheld from forward use because of fragmentation hazard and excessive weight, while portable oxygen concentrators fail without continuous electrical power or controlled environmental conditions. A novel approach is therefore required to enable safe, reliable oxygen therapy at or near the point of injury - the setting in which many preventable battlefield deaths occur.

2. PROBLEM STATEMENT / CAPABILITY GAP

Combat casualty care in large-scale combat operations (LSCO) and austere environments requires the ability to deliver oxygen therapy at or near the point of injury. Existing oxygen systems - including pressurized cylinders, portable therapeutic liquid oxygen (PTLOX) systems, and portable oxygen concentrators (POCs) - are constrained by hazard potential, power dependency, and logistical burden. These limitations impede delivery of lifesaving oxygen during the most time-sensitive phase of trauma care.

Pressurized oxygen cylinders, although capable of supplying high-purity oxygen, present fragmentation and blast hazards when exposed to small-arms fire, indirect fire, or thermal compromise. Owing to these risks, cylinders are frequently withheld from forward, dismounted, and high-threat operations, effectively denying oxygen therapy at the point where casualties are most likely to occur. Their mass and volume further restrict carriage on casualty evacuation (CASEVAC) platforms and within medical aid bags, forcing triage of oxygen resources and increasing the probability of hypoxic deterioration.

Portable oxygen concentrators require continuous electrical power and a clean, consistent air supply; performance degrades under conditions of altitude, high ambient temperature, humidity, dust load, and elevated respiratory demand (Allison & Dyer, 2018). In LSCO environments characterized by electromagnetic interference, degraded power availability, and dispersed formations, concentrators cannot be relied upon as the primary source of therapeutic oxygen.

During prolonged casualty care (PCC), extended evacuation timelines, terrain restrictions, airspace denial, and contested logistics compound the risk of hypoxia. Future forces will operate with greater dispersion, mobility, and airfield vulnerability, further reducing access to conventional oxygen systems. This combination of factors constitutes a persistent and well-documented capability gap: frontline medical personnel cannot safely or reliably deliver oxygen during the period when it is most critical to survival.

The absence of early oxygen administration contributes directly to preventable mortality from airway compromise, traumatic brain injury, and respiratory failure. Moreover, this capability gap is not confined to military operations. Civilian emergency medical services, firefighting units, search-and-rescue teams, and law-enforcement agencies encounter analogous constraints during disasters, structural collapses, and hazardous-materials incidents, where electrical power is unreliable and compressed-gas systems are present unacceptable safety risks. The deficiency in forward oxygen capability is therefore universal across both military and civilian emergency-response domains.

Problem Statement

The U.S. military and civilian responders lack a safe, portable, power-independent oxygen system capable of supporting respiratory interventions in the most dangerous and logistically constrained operating environments. This shortfall reduces survivability, limits TCCC execution, and exposes casualties to avoidable hypoxic injury.

3. OPERATIONAL IMPACT OF THE GAP

Hypoxia accelerates the transition from survivable injury to irreversible physiological collapse. In trauma casualties, inadequate oxygen delivery precipitates progressive shock, metabolic acidosis, coagulopathy, and secondary brain injury. Combat casualty data from Iraq and Afghanistan identify oxygen failure - resulting from hemorrhage or airway compromise - as a major contributor to potentially survivable deaths (Eastridge et al., 2012). Insufficient tissue oxygenation and systemic hypoperfusion promote trauma-induced coagulopathy and metabolic failure, amplifying mortality risk (Brohi et al., 2003). Even brief episodes of prehospital hypoxia have been shown to triple mortality in patients with traumatic brain injury and worsen neurologic outcomes (Spaite et al., 2017).

Across U.S. combat operations from 2001 to 2011, approximately 24.3 percent of battlefield deaths were assessed as potentially survivable, and 87 percent of those fatalities occurred before arrival at a medical treatment facility. Airway and respiratory compromise represented the second-leading cause of potentially survivable death - conditions in which early oxygenation remains a critical intervention (Eastridge et al., 2012).

Additionally, Damage Control Resuscitation (DCR) doctrine emphasizes the prevention of hypoxia, shock progression, and coagulopathy as core components of early trauma care. A central objective of DCR is the mitigation of oxygen debt—the cumulative physiological deficit that develops when oxygen delivery fails to meet metabolic demand. Even brief periods of uncorrected hypoxia accelerate metabolic acidosis, inflammatory cascade activation, and coagulopathy, driving the casualty toward irreversible physiological collapse. As multiple authors note, “it is not possible to incur a significant oxygen debt and suffer no consequences for lack of timely repayment,” underscoring that oxygen is not merely supportive but time-critical (Damage Control Resuscitation, 2023; Cap et al., 2015). In prolonged prehospital and contested evacuation environments, the inability to administer oxygen early directly contradicts DCR objectives and increases preventable mortality.

Future large-scale combat operations are expected to exacerbate these risks through delayed evacuation and prolonged exposure on the battlefield. Contributing operational factors include:

- Contested air corridors limiting medical evacuation (MEDEVAC) launch approval
- Extended ground evacuation distances resulting from distributed formations
- Degraded or denied communications impeding evacuation requests
- Urban density and rubble restricting vehicle access
- Chemical, biological, radiological, and nuclear (CBRN) threat environments necessitating protective posture

These conditions keep casualties on the X - the center of the fight - longer, forcing medics to manage critical injuries under continued threat and without immediate evacuation. Without reliable oxygen access, medics must manage casualties with respiratory distress, TBI, or thoracic injury under conditions that preclude prevention of hypoxic neurological decline.

The operational implications extend beyond the military. In joint civilian–military operations, including humanitarian assistance and disaster response, infrastructure degradation - such as power loss, debris entrapment, or toxic atmospheres - similarly restricts access to conventional oxygen systems. First responders, firefighters, and emergency medical services face the same functional gap, with corresponding adverse effects on patient survivability and neurological outcomes.

Operational Impact Statement


Extended hypoxia resulting from oxygen delivery shortfalls during prehospital care increases preventable mortality and disability, reduces combat power availability, and degrades mission readiness. Closing this gap is essential to achieving survivability gains in LSCO and high-risk domestic response operations.

4. OPERATIONAL & PHYSIOLOGICAL CONSEQUENCES OF HYPOXIA

Hypoxic deterioration occurs rapidly following traumatic injury. When respiration is impaired, inadequate oxygen delivery initiates a cascade leading to cellular dysfunction, secondary brain injury, multi-organ failure, and ultimately death. In traumatic brain injury (TBI), hypoxia and hypotension are the two most significant predictors of mortality and long-term neurological deficit (Spaite et al., 2017; Dutton et al., 2002). Even brief periods of oxygen deprivation substantially increase the likelihood of irreversible neuronal damage (Yan et al., 2014; Seo et al., 2019).

During prehospital combat care, hypoxia may occur as a result of airway obstruction, impaired ventilation, thoracic trauma, blast lung, smoke or toxic inhalation, or loss of consciousness following blunt or penetrating trauma. Without supplemental oxygen, casualties deteriorate along a predictable trajectory driven by shock, catecholamine depletion, and impaired perfusion (Eastridge et al., 2012; Kotwal et al., 2016). Timely interventions within the first hour post-injury - the so-called Golden Hour - are strongly associated with improved survival and neurological preservation (Howard et al., 2018; Mabry & Kotwal, 2019). This progression is illustrated in Figure 1.

Figure 1. Hypoxic Deterioration Sequence in Combat Trauma



Oxygen therapy is therefore not a supportive treatment - it is a primary life-saving intervention. The ability to administer oxygen immediately after injury is essential to prevent secondary complications, avoid respiratory arrest, and preserve neurological function until evacuation or advanced care is available.

Operational Reality: The highest proportion of preventable deaths occur in areas where conventional oxygen systems cannot be safely or reliably deployed.

5. EVALUATION CRITERIA FOR AN EFFECTIVE SOLUTION

Any oxygen delivery system intended for LSCO, austere, or hazardous environments must meet stringent performance criteria to ensure both clinical effectiveness and operational suitability.

5.1 Clinical Performance Requirements

A suitable system must:

- Deliver ≥ 90 % oxygen purity suitable for airway and ventilatory support (USAISR, 2018; JTS CPG Airway Management, 2020).
- Maintain consistent flow ≥ 6 L/min for therapeutic application (ATP 4-02.55, 2019; JTS CPG Airway Management, 2020).
- Enable support for airway compromise, TBI, burns, and thoracic trauma (CoTCCC Guidelines, 2023).
- Integrate with standard nasal cannula, masks, and ventilators (MIL-STD-1472G, 2012; Zoll Impact 731 Operator's Manual).
- Be immediately deployable with minimal cognitive and motor load on medics (ATP 4-02.55, 2019; HSI Standards for Medical Equipment, 2017).

These factors ensure the system can support TCCC Guidelines and prevent secondary hypoxic brain injury during prolonged prehospital care (JTS CPG TBI Management, 2020).

5.2 Operational Suitability Requirements

To function where casualties occur, systems must be:

- Non-pressurized to eliminate explosive fragmentation hazard
- Power-independent to ensure performance in degraded grids or blackout conditions
- Man-portable with minimal weight and space burden
- Ruggedized to withstand the environmental extremes of LSCO
- Zero-maintenance with long shelf life and no refill dependencies
- Approved for global transport via ground, air, and maritime supply networks
- Intuitive enough to require no specialized training or sustainment

These parameters reflect doctrinally informed requirements drawn from TCCC, JTS guidelines, and Army Health System operational expectations.

5.3 Survivability and Readiness Considerations

- Must be safe aboard rotary- and fixed-wing aircraft, tactical vehicles, naval platforms
- Must not contribute to logistics vulnerabilities (power, batteries, hazardous materials)
- Must remain functional aftershock, vibration, temperature swings, or contaminants
- Must align with future force concepts emphasizing dispersion, mobility, and contested logistics

6. SOLUTION OVERVIEW: RUGGED OXYGEN GENERATOR (ROG)

The ROG provides a solid-state, non-pressurized oxygen capability designed specifically for austere, high-risk, and power-denied environments. Unlike compressed gas or concentrator systems, the ROG™ utilizes a chemical oxygen generation reaction of sodium chlorate (2NaClO_3) to deliver a minimum of 96% oxygen at ≥ 6 L/min for a minimum of 15 minutes, a total of 90 liters of O₂ capacity, with no external power, batteries, or refills required.

6.1 System Design & Components

The system consists of:

- A sealed oxygen-generating chemical canister
- A simple twist-activation mechanism
- Standard oxygen outlet port
- Built-in Cool Touch thermal mitigation layer
- Ruggedized outer casing optimized for field handling

The system is entirely self-contained and maintenance-free throughout its operational life. It requires three actions to deploy:

Connect → Twist → Administer - enabling treatment in seconds with low cognitive demand.

6.2 Performance Parameters

- Purity: Minimum of 96% medical-grade oxygen
- Flow rate: ≥ 6 liters per minute
- Duration: ≥ 15 minutes
- Weight: ~2.45 kg (man-portable)
- Shelf life: 1 year minimum
- Deployment orientation: 360° usable position
- Interfaces with masks, cannula, and ventilators with standard low-flow bypass fittings


The system can be employed as:

- Primary oxygen capability for dismounted or squad-level care
- Supplemental oxygen for ventilated casualties (including during MEDEVAC)
- Emergency oxygen during power failure or equipment compromise

6.3 Operational Advantages

- Zero explosive hazard: safe in direct or indirect fire environments per DEF STAN 00-35 ballistics testing
- No power reliance: functions during blackout and silent operations
- Low signature: no emitted noise, or light detectable to the enemy
- Interoperability: compatible with current U.S. medical kits and platforms
- Extreme environment resilience: cold, heat, dust, altitude, humidity
- No roll design - allows the ROG to be set down during use
- No fire accelerators as any part of the system
- Low cognitive load: No adjustment valves, flow regulation tools, or power sequencing required

Figure 2. ROG Operational Advantages



6.4 Gas Conditioning & Thermal Moderation

The oxygen-generation process is exothermic; however, the ROG incorporates passive thermal moderation and gas-conditioning structures that ensure oxygen exiting the device is cooled, filtered, and safe for direct patient delivery at ambient breathing temperature.

Key characteristics:

- Thermal output is contained and moderated within the sealed oxygen canister and device housing.
- Oxygen is conditioned and stabilized prior to reaching the outlet.
- The exterior surface remains safe to handle throughout operation.
- No powered cooling or external equipment is required.
- This enables immediate clinical use without additional control hardware.

ROG provides immediate oxygen therapy at the tactical edge where traditional oxygen systems cannot be risked or sustained.

6.5 Role Within the Existing Medical Oxygen Ecosystem

The ROG is not intended to replace conventional oxygen systems such as compressed cylinders, PTLOX units, or powered concentrators. These systems remain essential for sustained, high-duration oxygen therapy at higher levels of care. Rather, the ROG™ addresses a capability gap in the earliest and most time-sensitive phase of casualty management, where traditional oxygen systems are either too hazardous, too heavy, or too power-dependent to employ.

ROG fills the operational requirement for:

- Immediate oxygen delivery at or near the point of injury
- Care in dismounted, remote, or confined environments
- Low-signature operations where power or pressurized systems pose risk
- Prolonged field care when evacuation is delayed or denied

In this role, the ROG complements existing medical oxygen infrastructure by extending oxygen availability forward - ensuring casualties receive critical oxygenation earlier, improving stabilization prior to handoff, evacuation, or escalation to higher medical roles.

7. TESTING & VALIDATION

ROG performance and environmental durability have been successfully verified through formal testing aligned with DoD medical system expectations.

7.1 MIL-STD-810G Environmental Qualification

The system has completed testing per MIL-STD-810G including:

Test	Test Standard
Vibration	MIL-STD-810G/H, Method 514.8, Procedure I
High Temperature	Tailored from MIL-STD-810G, Method 501.5, Procedure III
High Humidity & High Temperature	Tailored from MIL-STD-810G, Method 507.5
Low Temperature	Tailored from MIL-STD-810G, Method 502.5, Procedure II
Altitude	MIL-STD-810G, Method 500.5, Procedure II
Blown Sand	MIL-STD-810G, Method 510.6, Procedure II
Blown Rain	MIL-STD-810G, Method 500.6, Procedure I
Operational Drop Test	MIL-STD-810G, Method 516.6, Procedure IV
Mechanical Shock	MIL-STD-810G, Method 516.7, Procedure I

Table 1. Summary of MIL-STD-810G Tests Performed

ROG maintained consistent oxygen purity and flow throughout these environmental stresses. This testing ensures the ROG operates reliably in real battlefield heat, cold, dust, vibration, and impact conditions according to military standards.

7.2 FDA and Global Transport Certification

- FDA 510(k) cleared medical device
- Classified as UN3356 - permitted for global logistics via:
 - Commercial & military airlift
 - Ground and maritime transport
 - Storage aboard aircraft, ships, tactical vehicles

This eliminates hazardous materials restrictions that reduce availability of pressurized oxygen systems and verify the oxygen output is clinically safe and approved for patient treatment, not experimental use.

7.3 Ventilator Compatibility Testing

The ROG was validated with the Zoll Impact 731 Ventilator under simulated patient settings. This confirms stable oxygen delivery when used with standard ventilators, supporting seamless care during CASEVAC. Results demonstrated:

- Stable oxygen delivery of ≥ 6 L/min
- No oxygen purity degradation
- Ventilator alarms triggered only for incomplete exhalation, not oxygen failure
- No mechanical or safety conflicts documented

7.4 Quality Assurance & Production Standards

- Manufactured under a validated, ISO-compliant QA framework
- LOT traceability and sample testing ensure consistent output and reliability

7.5 Legacy System Validation and Proven Technology Base

The Rugged Oxygen Generator (ROG) builds on Molecular Products' decades-long expertise in chemical oxygen generation; the same core technology used in multiple fielded and certified systems across defense and industrial sectors.

- CAN-33 Chemical Oxygen Generator - Produces approximately 3,300 litres of breathable oxygen for submarine and mine-refuge applications. The CAN-33 has been integrated into U.S. and NATO naval platforms, demonstrating consistent oxygen purity and safety in enclosed environments.
- Multi-Purpose Oxygen Generator (MPOG) and Electric MPOG (eMPOG) - These systems provide up to 2,600 litres of oxygen per unit for medical and life-support contingencies. Their sustained success under high-humidity, vibration, and temperature extremes validates the stability and controllability of the underlying chemical oxygen formulation.

Proven Field Performance - These legacy systems have logged tens of thousands of operational hours across maritime and industrial platforms with zero catastrophic failures.

This proven heritage provides a verified technical foundation for the ROG. By adapting the same stabilized oxygen-generation chemistry into a lightweight (<3 kg), power-independent design, the ROG transitions a trusted maritime oxygen technology into a deployable capability for land-based medical and combat operations.

8. COMPARATIVE ANALYSIS

To evaluate the operational effectiveness of the ROG, its performance must be compared against other field oxygen delivery systems currently available to military and civilian responders. The most widely utilized systems - compressed gas cylinders and portable oxygen concentrators (POCs)

each carry limitations that restrict safe deployment into the operational environments where oxygen therapy is most needed.

8.1 Limitations of Compressed Oxygen Cylinders

- Fragmentation and blast hazard: catastrophic risk under weapons fire
- Weight and bulk: limits portability and storage aboard CASEVAC platforms
- Requires refill supply chain, adding logistics burden in LSCO
- Often restricted from dismounted operations due to blast hazard risk
- Potential to become inadvertent missile hazard in vehicle roll over or aircraft mishap

Even where available, cylinders may be intentionally withheld from high-threat areas, preventing use at or near the point of injury.

8.2 Limitations of Portable Oxygen Concentrators

- Require continuous power (batteries or generator-based systems)
- Degraded oxygen purity and flow at altitude or in thermal extremes
- Dust and humidity ingress reduce function in austere terrain

- Noisy mechanical operation can compromise signature management
- Lack ruggedization expected for LSCO environments

In degraded operational conditions - the exact scenarios where hypoxia risk increases - POCs become unreliable or non-functional.

8.3 System-Level Comparative Table

Capability Feature	ROG	Compressed Cylinders	POCs
Non-pressurized (no blast hazard)			
Requires electrical supply			
High oxygen purity maintained in extreme environments			
Man portable for dismounted use			
Shelf-life, no refills or charging			
Silent / no electro-optical signature			
Shelf-stable, low sustainment burden			
Safe for aircraft and tactical vehicles			
Rapid deployment under fire			

Possible but commonly constrained by hazard, power, or sustainment limitations

Table 2. Capability Comparison of Oxygen Systems

Conclusion

ROG is the only oxygen system safe and reliable enough to support dismounted, contested logistics, and high-threat tactical environments.

9. CROSS-DOMAIN MEDICAL APPLICATIONS

While developed to meet combat medicine requirements, the ROG solves a multi-agency capability gap shared by national and local emergency responders.

9.1 Military Use Cases

- Squad-level and point-of-injury care where cylinders are not authorized
- Prolonged Casualty Care (PCC) during delayed evacuation
- Airway and ventilatory support during ground and air CASEVAC
- Tactical operations inside confined spaces, subterranean terrain, or high-rise structures
- CBRNE / toxic exposure response where ambient air cannot be used

ROG expands oxygen availability deep into the fight, directly supporting TCCC execution.

9.2 Military Use Cases

ROG provides essential capability for:

- EMS operating in remote or disaster zones with no power access
- Firefighters entering oxygen-depleted or toxic environments
- Search & rescue during building collapse or industrial accidents
- Law enforcement in barricade, overdose, or hazmat situations
- Rural healthcare, where transport distances prolong hypoxic risk

In these use cases, lives are lost when oxygen cannot be administered early - a gap ROG uniquely fills.

9.3 Joint Operations and Interoperability

Because oxygen delivery standards are universal, ROG requires:

- No specialized training
- No device-specific interfaces
- Minimal kitting changes
- Pictogram instructions and tamper proof cover latch

This interoperability increases readiness across:

- DoD medical forces
- Department of Homeland Security(DHS) components
- Federal/state/local first responders

Bottom Line:

ROG closes a shared operational and medical gap across both national defense and civil emergency care.

10. ACQUISITION AND SUSTAINMENT PATHWAY

The ROG is a commercial off-the-shelf (COTS) medical device that requires no new infrastructure, minimal training, and limited sustainment investment. Its adoption can be rapidly executed through existing U.S. military and federal procurement channels.

10.1 Acquisition Readiness

ROG meets the critical acquisition prerequisites for defense medical integration:

- Fully COTS available and in active production
- FDA 510(k) cleared for patient use
- UN3356 certified for global transport
- No hazardous materials licensing required for storage, transport, or disposal
- MIL-STD-810G validated for environmental ruggedness

These factors reduce technical risk and allow rapid transition from evaluation to fielding.

10.2 Acquisition Readiness

ROG can be procured through:

- DLA Medical supply chains(pending NSN assignment)
- GSA schedules for interagency access (pending NSN assignment)
- Multi-year sustainment agreements for unit readiness
- Rapid acquisition pipelines for urgent operational needs

No modifications to current medical equipment lists, evacuation kits, or clinical workflows are required.

10.3 Logistics and Sustainment Advantages

- No batteries or power accessories to replace
- No refills required throughout shelf-life
- No periodic maintenance or calibration
- Shelf-stable for one year minimum
- Reduced HQDA Class VIII burden during LSCO
- Eliminates bulk supply distribution of compressed oxygen to the edge

This logistics simplicity directly enhances operational reach and medical endurance under contested supply conditions.

10.4 Training and Integration

- Intuitive twist to-activate deployment
- Compatible with existing ventilation and airway equipment
- Minimal user training (~3 minutes)
- Easily integrated into TCCC and PCC curricula

Recommendation Summary:

The ROG provides an immediate, zero-power, non-hazardous oxygen capability suitable for the most dangerous and logistically constrained operating environments. Its adoption will reduce preventable deaths, increase force survivability, and improve readiness across defense and civilian emergency response.

11. REFERENCES

- Allison MG, Dyer RS. Wilderness environmental medicine considerations for prolonged field care. *Wilderness Environ Med.* 2018.
- ATP 4-02.55. Army Health System Support to Maneuver Forces. Department of the Army; 2019.
- Brohi K, Singh J, Heron M, Coats T. Acute coagulopathy of trauma: Hypoperfusion induces systemic anticoagulation. *J Trauma.* 2003;54(6):1127–1130.
- Cap AP, et al. “Blood failure”: Time to view blood as an organ—How oxygen debt contributes to blood failure and its implications for remote damage control resuscitation. *Transfusion.* 2015.
- Committee on Tactical Combat Casualty Care (CoTCCC). TCCC Guidelines. Defense Health Agency; 2023. Damage Control Resuscitation. Defense Health Agency; 29 Aug 2023.
- Dutton RP, Stansbury LG, LeBlanc LH, et al. Hypotension and hypoxia in major trauma: Causes of preventable death in the first hour. *J Trauma.* 2002;53(5):S183–S189.
- Eastridge BJ, Mabry RL, Seguin P, et al. Death on the battlefield (2001–2011): Implications for the future of combat casualty care. *J Trauma Acute Care Surg.* 2012;73(6 Suppl 5):S431–S437.
- Human Systems Integration Standards for Medical Equipment. Department of Defense; 2017.
- Joint Trauma System (JTS). Clinical Practice Guideline: Prehospital Airway Management. Defense Health Agency; 2020.
- Joint Trauma System (JTS). Clinical Practice Guideline: Traumatic Brain Injury Management. Defense Health Agency; 2020.
- MIL-STD-1472G. Department of Defense Design Criteria Standard: Human Engineering. Department of Defense; 2012.
- OpenCriticalCare.org. Oxygen delivery fundamentals; 2024.
- Seo DE, Shin SD, Song KJ, Ro YS, Hong KJ, Park JH. Effect of hypoxia on mortality and disability in traumatic brain injury according to shock status: A cross-sectional analysis. *Am J Emerg Med.* 2019;37(9):1709–1715.
- Spaite DW, Hu C, Bobrow BJ, et al. The effect of combined out-of-hospital hypotension and hypoxia on mortality in major traumatic brain injury. *Ann Emerg Med.* 2017;69(1):62–72.
- U.S. Army Institute of Surgical Research (USAISR). Oxygen Delivery Systems in Prolonged Field Care. 2018.
- Yan EB, Satgunaseelan L, Paul E, et al. Post-traumatic hypoxia is associated with prolonged cerebral cytokine production, higher serum biomarker levels, and poor outcome in severe TBI. *J Neurotrauma.* 2014;31(7):618–629.
- Zoll Medical Corporation. Impact 731 Ventilator Operator’s Manual. 2019.